



Real-Time Actionable Data After Each UVC Room Disinfection

Introducing DoseEZ™ by ConfluenceAI - a patented technology designed to improve infection control in healthcare settings. These advanced UVC sensors provide real-time data after each disinfection, helping to combat Hospital Acquired Infections (HAI). With cloud storage of disinfection treatment data, DoseEZ is a powerful tool in the fight against HAI.


One DoseEZ

**USE WITH ANY UVC UNIT
UNLIMITED TREATMENTS
ENTIRE HOSPITAL
JUST \$199/MONTH**


ANNUAL SUBSCRIPTION LAST MONTH FREE (\$182/Mo)



U.S. Patent No. 12,491,282-B2
US and Int'l Patents Pending



Effect of UVC Device Placement on Angle
Compare the Angles to UVC Reduction in Fig. 2
Fig. 1



Horizontal Surface Patient Bed or OR Table
85 Degree Angle = 9% UVC Energy
Fig. 2

LAMBERT'S COSINE LAW
UVC LIGHT EFFECTIVENESS
DRAMATICALLY LOWERS AS THE
ANGLE INCREASES

If your UVC Device reports 100mJ/cm² UVC Dose when shining directly onto a surface (Fig 2. 100%), then at a 60 degree angle, your OR table or patient bed is only receiving 50% (50mJ/cm²)!

DoseEZ compensates for angled UVC for the highest accuracy. Ignore your UVC Device reports and use DoseEZ's comprehensive data.


Infection Prevention Managers

Customize The Cloud Dashboards To Your Needs

OR12		OR12 Dose	220 mJ/cm ² Last Dose	OR12 Threshold — 4 days ago	Max: 200	OR12 Results — 4 days ago	1.92 min last_or12_duration
OR16		OR16 Dose	126 mJ/cm ² Last Dose	OR16 Threshold — 4 days ago	Within Range	OR16 Results — 4 days ago	2.18 min last_or16_duration
OR17		OR17 Dose	9 mJ/cm ² Last Dose	OR17 Threshold — 4 days ago	Within Range	OR17 Results — 4 days ago	01/04 13:58 EST last_or17_treatment
OR6		OR6 Dose	225 mJ/cm ² Last Dose	OR6 Threshold — 4 days ago	Max: 200	OR6 Results — 4 days ago	01/09 15:06 EST last_or6_treatment

Easy to Understand Color Coded Results

- This table includes:
- Days Since Last Treated
- Dose Amount:
 - Yellow <10% Overage
 - Green Within Limits
 - Rose Under Limit or >10% Over
- Threshold Shows Min & Max Values or if Within Range
- Results Shows Treatment Duration and Date/Time
- Custom Reports Available



Choose Various Data Logs, Charts and Widgets

This table includes many valuable datapoints, such as Building, Room, Dose, Duration, Alert Type, Maximum UVC Treatment Level, Minimum & Maximum UVC Limits & more.

Goodbye Dots, Hello IoT Sensors!

DoseEZ replaces the expensive color-changing paper dots with a digital solution that measures precise dose levels, detects over/under-dosing, and permanently documents historical records.



Validate Treatment Dose

Patented DoseEZ uses QR code identifiers to Cloud-Store UVC doses for EACH sensor location in EVERY room. Add QR codes to each UVC device. The Cloud will use historical performance data, and alert staff when new lamps are needed.



Immediately Actionable Cloud Data

Dose treatment levels sent to the Cloud will immediately alert staff members, if the UVC dosage is below or above your set limits.

Nothing To Purchase

DoseEZ is an affordable monthly or annual subscription and not a capital equipment purchase. The subscription includes hardware maintenance and software updates.

Specifications

- Wireless Connectivity: Bluetooth® and WiFi. Optional LoRaWAN Eliminates WiFi Need.
- No HIPAA Issues. Use Hospital Guest WiFi network
- Data Storage: Cloud Enterprise Grade Management Environment
- Phone Apps: IOS and Android
- Power: 12VDC. Rechargeable Battery for Mobile Use
- NIST Traceable 254nm UVC Sensor